# CHAPTER III EXPERIMENTAL

## 3.1 Materials

- 1. Sodium hydroxide (NaOH, Labscan Asia Co.)
- 2. Sulfuric acid (H<sub>2</sub>SO<sub>4</sub>, Merck Co., Germany)
- 3. Hydrochloric acid (HCl, Labscan Asia Co.)
- 4. High purity standards: Glucose (Sigma Aldrich Chemical Co. Inc., USA)
- 5. High purity standards: Xylose (Sigma Aldrich Chemical Co. Inc., USA)
- 6. High purity standards: Arabinose (Sigma Aldrich Chemical Co. Inc., USA)
- 7. Calcium carbonate (CaCO<sub>3</sub>, Sigma Aldrich Chemical Co. Inc., USA)
- 8. Sodium sulfite (Na<sub>2</sub>SO<sub>3</sub>, Labscan Asia Co.)
- 9. Cellulase from *Trichoderma reesei ATCC 26921* (Sigma Aldrich Chemical Co. Inc., USA)
- 10. Yeast extract (Himedia Laboratories, India)
- 11. Malt extract (Himedia Laboratories, India)
- 12. Peptone (Himedia Laboratories, India)
- 13. Baker's yeast, *Saccharomyces cerevisiae* (Thailand Institute of Scientific and Technological Research (TISTR), Thailand)

# 3.2 Equipment

- 1. Ethos Sel: Microwave solvent extraction lab station (Milestone)
- Shimadzu Corp.: High Performance Liquid Chromatography (HPLC) with a refractive index detector(RID-10A, Shimadzu Corp., Kyoto, Japan) using a chromatography column Aminex-HPX 87H column (300 mm x78 mm, Bio-Rad Lab, USA) under 20 μl injection volume, 0.005 M sulfuric acid (HPLC grade) as

mobile phase, 0.6 ml/min flow rate, 65 °C column temperature, and 20 min run time condition

- 3. Scanning Electron Microscope (SEM) (Hitachi/S-4800)
- 4. Agilent Technologies: Gas Chromatography (GC)
- 5. Microscope
- 6. Analytical balance
- 7. Oven
- 8. Water bath
- 9. Autoclave
- 10. Rotary evaporator
- 11. Shaking Incubator

## 3.3 Methodology

## 3.3.1 Mission Grass preparation

Mission grass (*Pennisetum polystachion*) was harvested in March 2013 at Tak Province. The grass (stems and leaves) was left to dry under sunlight for 2 weeks. Then, it was milled into small particles (diameter  $\sim 300~\mu m$ ) by a milling factory. The milled mission grass was stored in a dry container at room temperature.

#### 3.3.2 Pretreatment of Mission Grass

Mission grass was pretreated with the grass optimization method according to Boonmanumsin et al. and Tatijarern et al. The milled grass was mixed with 3% w/v NaOH using liquid-to-solid ratio, 15 mL of NaOH: 1 g of grass (Boonmanumsin *et al.*, 2012, Tatijarern *et al.*, 2013). The mixture was stirred thoroughly until it was homogenous. Then, it was subjected to microwave treatment (300 W) under the

temperature of 120 °C for 10 minutes. After the pretreatment, the mixture was washed until the pH was neutral. The solid residue was oven dry for 24 hours.

#### 3.3.3 Hydrolysis of Mission Grass

The dry solid from the base pretreatment was mixed with 1% v/v H<sub>2</sub>SO<sub>4</sub> with liquid-to-solid ratio 15 mL of H<sub>2</sub>SO<sub>4</sub>: 1 g of pretreated grass. The mixture was stirred thoroughly, and then subjected to microwave treatment (300 W) under the temperature of 200 °C for 5 minutes (Tatijarern *et al.*, 2013). The mixture was left to cool at room temperature. Cellulase from *Trichoderma reesei ATCC 26921* (160 μL/1 g of grass) was added into the liquid hydrolyzate and incubated for 60 hours at 50°C (Qureshi *et al.*, 2008).

## 3.3.4 Detoxification of Mission Grass Hydrolyzate

Overliming method was used to detoxify the liquid hydrolyzate. The hydrolyzed solution was treated with Ca(OH)<sub>2</sub> to increase the pH to pH 8-12, then sodium sulfite (1 g/L) was added into the solution. The overlimed solution was heated for 30 minutes at 90°C, and the pH was adjusted to 6 using 95% H<sub>2</sub>SO<sub>4</sub> at room temperature. The condensate in the mixture was filtered. The liquid was stored for fermentation.

## 3.3.5 Fermentation of Mission Grass Hydrolyzate

# Preparation of cell culture

To prepare the yeast medium (YM) for *Saccharomyces cerevisiae*, yeast extract (3 g/L), malt extract (5 g/L), peptone (5 g/L), and glucose (20 g/L) were mixed and stirred until homogenous (Huang *et al.*, 2009). The medium was autoclaved at 121 °C for 15 minutes and cooled down to room temperature.

Agar plate medium was prepared by mixing glucose (10 g/L), yeast extract (3 g/L), malt extract (5 g/L), peptone (5 g/L), agar (20 g/L), and stirred at 100°C until homogeneous. The solution medium was autoclaved at 121°C and then poured into a petri dish under a biosafety cabinet. The agar in the plate was left to cool at room temperature for one day.

A vial of yeast strain obtained from TISTR was resuspended in a small amount of YM medium. A sanitized inoculation loop was dipped into the solution and streaked onto the prepared agar plate. The plate was incubated at 30°C for 24 hours.

A Saccharomyces cerevisiae. colony from the agar plate was chosen and transferred to the prepared YM medium. The inoculated medium was incubated in a shaker at 30°C for 24 hours. The agitation of the incubator was set to 100 rpm.

#### Inoculation of Saccharomyces cerevisiae

Yeast extract (3 g/L), malt extract (3 g/L), and peptone (5 g/L) were added into the overlimed grass hydrolyzate. No additional sugar was added. The mixture was then autoclaved at 121°C for 15 minutes and allowed to cool down to room temperature. Saccharomyces cerevisiae (10% v/v) was added into the autoclaved hydrolyzate. The solution was incubated in the shaker at 30°C for 96 hours. Sugar utilization and ethanol production was detected using HPLC and GC, respectively, every 24 hours.